

Natural detergent of vegetable origin



KOKONAT 2.0 - Endo Natural Cleaner

- Natural
- Effective
- Easy rinsing
- Contains no alcohol

DESCRIPTION

KOKONAT 2.0 is a natural cleaner based on anionic surfactants of plant origin obtained through biotechnology. The lowering of the surface tension of the liquid and the formation of micelles effectively remove dirt from surfaces. KOKONAT 2.0 can be used on any surface and material. It rinses off easily with a little water and contains no alcohol.

COMPOSITION

<i>Ingredient</i>	<i>% w/w</i>
Water	82.6
Anionic surfactants of vegetable origin	1.4
Carboxylates	16.0

APPLICATIONS

Manual and automatic cleaning of medical devices and their accessories (endoscopes, surgical instruments, etc.). The product can also be used in ultrasonic devices.

MODALITIES OF USE

Washing by hand or in ultrasonic tanks

Dilute 1:500 (2ml per litre of water). Immerse the device in the solution for 5 minutes; if there are persistent and/or long-dried organic residues, double the concentration and immersion time.

Rinse with water.

The product can be diluted in tap water.

CONSERVATION

Keep product in original containers. Store in a cool, dry place, at a temperature between 5 and 30 °C, away from any source of heat and direct exposure to sunlight.

Keep containers tightly closed. Ensure adequate room ventilation.

The product retains its effectiveness for up to 24 months from the date of production.

Natural detergent of vegetable origin

SAFETY AND ENVIRONMENTAL IMPACT

Avoid contact with eyes and skin.

Non-flammable product.

See safety data sheet for further information.

PACKAGING

Order Code	Model
DT-KKN2-5000	5 litre canister
-	Pack of 4 x 5-litre canisters

DISPOSAL OF THE DEVICE

Subject to local regulations. Do not dispose of container in the environment. The container cannot be reused.

DIRECTIVES AND COMPLIANCE

This product complies with the requirements of Regulation (EU) 2017/745 on medical devices.

Classification according to Regulation (EU) 2017/745 Article 51 - Annex VIII rule 1

Class I

CND: Z12029085
RDM: 2432134/R
Basic UDI-DI: 805053900DT-KKN266
UFI: FR10-90JY-FOOD-2EWX